


Kidsoft[®] Laparoscopic Pyeloplasty Stent Kit (Ureteral Stent & Set)

Intended Use:

Implanted for drainage of urine from renal pelvis to urinary bladder.
(After pyeloplasty procedure)



Size (in inch)

2.8, 3, 3.5, 4, 5 & 6

Length (in cm)

12, 14, 15, 16, 18 & 20

Set Includes

Introducer with Needle, Guidewire, Stent and Pusher

Key Features and Benefits:

- Ergonomically designed for stent placement during laparoscopic pyeloplasty
- The needle can be percutaneously introduced and placed across the anastomosis after removing the sharp stylet
- The stent can be guided into the ureter with ease, the proximal end can be brought out as a nephrostomy if introduced through the Kidney or Pelvis
- The stent can be retained for upto 30 days