

Malecot Catheter

Intended Use:

Malecot Catheter is intended to be inserted percutaneously into the renal pelvis to perform external drainage of urine



Package Includes

Malecot Catheter, metal or plastic Straightner and urine bag connector

Size (in Fr)

8, 10, 12, 14, 16, 18, 20 & 22

Length (in cm)

20, 25, 30 & 35

Key Features and Benefits:

- The Malecot flower design is suitable for drainage of thick viscous fluids
- Integral flower design matches with the catheter size when closed and cause less trauma during insertion