


Mono J Catheter

Intended Use:

Used for temporary drainage and irrigation of urine from the kidney (Pigtail ureteral catheter).



Types

Open & Closed End

Size (in Fr)

4.7, 5, 6, 7 & 8

Length (in cm)

70

Also available in set, which Includes
Mono J stent and guidewire

Key Features and Benefits:

- Pigtail loop helps retention in renal pelvis
- Holes in pigtail loop helps in better drainage
- Used to inject contrast media
- Used for short term external drainage to avoid stenting