

## N Quadra® (Nitinol Stone Basket With Tip)

### Intended Use:

Stone Basket is intended to be used in an endoscopic procedure for removal of urinary stones or foreign objects from the Ureter or the Kidney



**Size (in inch)**  
2.5, 3, 4 & 5

**Length (in cm)**  
70, 90, 110 & 120

### Key Features and Benefits:

- Nitinol wire construction with high strength and memory helps the basket to retain its shape during extreme torsion
- In case of stone impregnation in the basket during the procedure, ergonomic handle permits full disassembly
- Atraumatic tip