

## Nephrostomy Balloon Dilator

### Intended Use:

Used for dilation of percutaneous tract for nephrostomy procedure



### Product Specifications

- Shaft size: 5 Fr
- Balloon expansion Diameter : 10 mm
- Balloon length : 12 cm
- Catheter length : 55cm
- RBP : 11 atm
- Use guidewire : 0.035 inch

### Key Features and Benefits:

- Controlled inflation pressure allows dilating force at a stricture or narrowing
- Radiopaque marker bands enhance visualization of balloon position under fluoroscopy