


**PFlex**  
Ureteral Stent  
**Double J Stent**

**Intended Use:**

Intended to be implanted for drainage of urine from renal pelvis to urinary bladder



**Types**

Both Ends Open, One End Closed, Both Ends Multiloop, One End Multiloop  
With or Without Hydrophilic coating, With or Without String

**Package Includes**  
Stent with pusher

**Also available in set, which Includes**  
Stent with pusher and guidewire

**Key Features and Benefits:**

- For indwelling up to 30 days
- Polyurethane based stent has uniform diameter and smooth surface reducing encrustation
- Optimal loop shape and high memory retention reduces the risk of migration
- Optimal radiopacity for better visualization

Type	Size (in Fr)	Length (in cm)
Standard	4.5, 4.7, 4.8, 5, 5.5, 6 6.5, 7 & 8	12, 14, 16, 18, 20, 22, 24 26, 28 & 30
Multiloop	4.7, 4.8, 5, 6, 7 & 8	18-24 & 24-30
Pediatric	2.8, 3, 3.5, 3.7, 3.8 & 4	6, 8, 10, 12, 14 & 16