


PCN Catheter With Locking Mechanism

Intended Use:

PCN Catheter is intended to be inserted percutaneously into the renal pelvis to perform external drainage of urine. The locking mechanism prevents accidental slippage of the catheter.



Package Includes

PCN Catheter with locking mechanism, metal or plastic straightener and urine bag connector

Size (in Fr)

6, 7, 8, 10, 12, 14 & 16

Length (in cm)

20, 22, 25, 30, 35 & 40

Key Features and Benefits:

- The pigtail tip configuration helps in the retention of the distal tip in the renal pelvis
- The locking mechanism prevents accidental slippage of the catheter
- Eye holes for better drainage
- Also available in Hydrophilic coating for easy insertion

Note: Once the catheter is deployed, pull the string, wind and lock to the hub of the catheter