


PCN Catheter

Intended Use:

PCN Catheter is intended to be inserted percutaneously into the renal pelvis to perform external drainage of urine.



Package Includes

PCN Catheter, metal or plastic straightener and urine bag connector

Size (in Fr)

6, 7, 8, 8.5, 9, 10, 12, 14, 16 & 18

Length (in cm)

20, 22, 25, 30, 35 & 40

Key Features and Benefits:

- The pigtail tip configuration helps in the retention of the distal tip in the renal pelvis
- Eye holes for better drainage
- Also available with Hydrophilic coating for easy insertion