

Prime Wire[®]

PTFE Nitinol Guidewire

Intended Use:

Used to gain access to desired anatomical location and to assist in placement of medical devices during diagnostic or interventional procedures



Types	Size (in inch)	Length (in cm)
Straight Tip	0.025, 0.032, 0.035 & 0.038	150

Key Features and Benefits:

- **All-in-one design:** Streamlines urology procedures
- **Hydrophilic floppy tip (5cm):** Facilitates access in tortuous anatomy
- **Teflon coating:** Enables low-friction navigation and easy device advancement
- **Single wire access:** Avoids wire exchanges; enables access and stenting with one wire
- **East access:** Reduces procedural time