

Relief[®] CIC

Hydrophilic Intermittent Catheter

Your Intermittent Self-use
Catheter for Ultimate **Comfort**
and **Convenience**



**NO
LUBRICANT GEL
REQUIRED**



With Protective Sheath
For safe
storage after use

Hydrophilic Coated
For easy insertion & removal,
without using Lubricating Gel

Comfortable

- Perfect balance of catheter flexibility and stiffness that prevents bending and promotes effortless insertion
- Blunt tapered tip and smooth eyelets minimizes trauma during insertion
- Lubrication rendered by hydrophilic coating that saves preparation time and cost

Convenient

- Ready-to-Use, Hydrophilic Catheter for a quick and easy insertion, after dipping in clean water
- Ergonomic gripper design bestows non-slip control during insertion and removal
- Multi-coloured funnels aid easy identification of sizes


Relief[®] CIC

Hydrophilic Intermittent Catheter

Key Features


Benefits

Highly lubricious,
sustained hydrophilic coating




Helps in easy insertion
and removal

Clear catheter tube
with balanced flexibility
and stiffness



Prevents bending
of tube with enhanced visibility

Smooth eyelets at
opposite orientations



Optimum drainage

Blunt tapered tip



Minimizes trauma
during insertion

Ergonomic gripper design



Better control during
insertion and removal

Flexible funnel



Provides control to
prevent urine spillage

BLUENEEM[®]
Clinicians Trust. Patients Wellbeing



Blue Neem Medical Devices Pvt Ltd
Plot Nos 270 & 271, Road No 5
Harohalli Industrial Area, Phase II
Kanakapura Taluk
Ramanagara Dist - 562112
Karnataka, India



+91 97399 72855



enquiry@blueneem.com



www.blueneem.com

Relief [®] - Hydrophilic Intermittent Catheter			
Male		Female	
Size (in Fr)	Length (in cm)	Size (in Fr)	Length (in cm)
10	35	10	15
12	35	12	15
14	35	14	15
16	35	16	15
18	35	18	15

360° Product
Experience

