

Single Step[®] Ureteral Dilator

Intended Use:

Used for one step dilation of the ureter through natural orifice prior to a diagnostic or interventional procedure



Size (in Fr)

8/6, 10/5, 10/6, 10/8, 12/5, 12/6, 12/7, 12/8
14/5, 14/6, 14/7, 14/8, 16/6 & 16/7

Length (in cm)

70

Key Features and Benefits:

- Long tapered tip avoids use of sequential dilation
- Smooth surface for atraumatic entry during dilation
- Also available with hydrophilic coating for ease of introduction
- Optimal radiopacity for better visualization