


Suprapubic Catheter

Intended Use:

Suprapubic Catheter is intended to be inserted percutaneously into the urinary bladder and perform external drainage of urine from the urinary bladder.



Package Includes

Suprapubic Catheter, metal straightner, metal trocar and urine bag connector

Size (in Fr)

8, 10, 12, 14, 16, 18 & 20

Length (in cm)

25 & 30

Key Features and Benefits:

- Integral flower design matches with the catheter size when closed and cause less trauma during insertion
- Trocar and catheter design to prevent buckling of catheter during the procedure