

Sure Puncture Needle (Initial Puncture Needle - Three Part)

Intended Use:

Used to make the initial percutaneous puncture to create access and to enable the introduction of a guidewire during urological and percutaneous interventional procedures.



| Type | Size (in G) | Length (in cm) |
|------------|-------------|----------------|
| Three part | 21/16 | 20 |

Key Features and Benefits:

- **Siliconized needle:** Reduces trauma for a smoother procedure
- **Echogenic tip:** Enhances visibility under ultrasound for precise guidance
- **1 cm marking:** Ensure accurate depth measurement
- **21 G bevel tip + 16 G cannula needle:** Provides optimal puncture balance for secure access
- **Inner 21 G needle:** Allows minimal site puncture and contrast injection
- **Coaxial needle:** Ensures stable and secured access
- **21 G puncture for guidewire insertion:** Facilitates smooth guidewire introduction