

## Ureteral Balloon Dilator

### Intended Use:

Used for dilation of narrow ureteral segments



### Product Specifications

- Shaft size: 5 Fr
- Balloon Diameter: 5 mm
- Balloon Length: 4 cm
- Catheter Length: 75 cm
- RBP: 14 atm
- Use Guidewire: 0.035inch

### Key Features and Benefits:

- Controlled inflation pressure allows dilating force at a stricture or narrowing
- Radiopaque marker bands enhance visualization of balloon position under fluoroscopy