


Urethral Dilator Set - S Curved

Intended Use:

Used for dilation of the Male Urethral strictures



Type
S Curved

Size (in Fr)
8-24
(8, 10, 12, 14, 16, 18, 20, 22 & 24)

Length (in cm)
40

Key Features and Benefits:

- Tapered end for atraumatic entry during dilation
- Smooth surface for ease of introduction
- Optimal radiopacity for better visualization
- To be used for sequential dilation
- Side hole for drainage
- S curve supports natural curvature of male urethra