

# Urethral Dilator Set - Straight

**Intended Use:**  
Used for dilation of the Urethra



<b>Type</b> Straight	<b>Size (in Fr)</b> 8-24 (8, 10, 12, 14, 16, 18, 20, 22 & 24)  8-30 (8, 10, 12, 14, 16, 18, 20, 22, 24, 26, 28 & 30)	<b>Length (in cm)</b> 40
<b>Key Features and Benefits:</b> <ul style="list-style-type: none"><li>• Tapered end for atraumatic entry during dilation</li><li>• Smooth surface for ease of introduction</li><li>• Optimal radiopacity for better visualization</li><li>• To be used for sequential dilation</li><li>• Side hole for drainage</li></ul>		